

# Evaluation of Nanospray Voltage and Spray Stability and their Impact on Chromatographic Peak Area

Gary Valaskovic<sup>1</sup>, Mike Lee<sup>2</sup>, Amanda Berg<sup>1</sup>

<sup>1</sup>New Objective, Inc., Woburn, MA <sup>2</sup>Milestone Development Services, Newtown, PA

## Introduction

State-of-the art liquid chromatography-tandem mass spectrometry (LC-MS/MS) analysis uses a constant electrospray (ESI) voltage for data acquisition. Modern qualitative and quantitative LC-MS/MS methods depend on highly efficient gradient elution chromatography. The changing chemical composition of mobile phase during gradient elution results in an inherent disconnect with single point ESI voltage optimization. Using a nanospray source equipped with a digitally controlled stage and software for precise and reproducible emitter positioning for data acquisition we investigate the relationship between spray stability and data quality. Repetitive on-column injections at different (fixed) target ESI voltage settings were executed for four separate analytes with a wide range of gradient elution times. Plotting the chromatographic peak area for selected ion currents yields an apparent compound dependant response curve in which a total maximum value is observed. Image capture enabled by the Digital PicoView software program reveals a direct correlation between the observed spray mode, spray stability and chromatographic quality.

## Methods & Materials

### Mass Spectrometer

- LTQ Linear Ion Trap (Thermo)
  - Full Scan MS: 150 -2000 Da
  - Spray Voltage: fixed per data file, variable across replicated injections
  - Analyte Specific Targeted MS Scans
    - MRFA 524 Da MH<sup>+</sup>: Mass Range 514.3–534.3 Da
    - Bradykinin 1-7 Fragment 379 Da MH<sub>2</sub><sup>+</sup>: Mass Range 374.4–384.4 Da
    - Bradykinin 531 Da MH<sub>2</sub><sup>+</sup>: Mass Range 521.0–541.0 Da
    - Neurotensin 837 Da MH<sub>2</sub><sup>+</sup>: Mass Range 827.2–847.2 Da
    - DPV-550 Digital PicoView nanospray source (New Objective, Inc.)

### Chromatography





- Agilent 1100 Capillary Pump
  - 20 µL Flow meter, Operated in MicroMode
  - Flow Rate: 1.0 µl/min
  - Mobile Phase A: 0.1% formic acid in water
  - Mobile Phase B: 0.1% formic acid in acetonitrile
  - Gradient: 30 minutes 2-50% B
- Column: PicoFrit column (360 µm OD x 75 µm ID x 15 µm tip) slurry packed to 10 cm with BioBasic C18 (5µm, 300 Å, C18, Thermo)
- HTC Pal Autosampler (Leap Technologies)
  - 6-port injection valve (VICI Valco Instruments, Inc.)
  - 1.0 µL loop

### Reagents







- MRFA: 500 fmol/µL, MW 523 Da (Sigma Aldrich)
- Bradykinin 1-7 fragment: 500 fmol/µL, MW 756.4 Da (Sigma Aldrich)
- Bradykinin: 500 fmol/µL, MW 1060.2 Da (Sigma Aldrich)
- Neurotensin: 500 fmol/µL, MW 1672.0 Da (Sigma Aldrich)
- 0.1% Formic acid in water (JT Baker)
- 0.1% Formic acid in acetonitrile (JT Baker)

## Conclusions

- Effect of voltage related spray stability on data quality across a gradient demonstrated
  - Gradient associated changes in spray stability observed
  - Spray voltage optimum according to maximum peak area observed
  - Increased noise and decreased peak area observed for poor spray stability
  - Ion signal stability appears to correlate with periodic (1-10 Hz) events associated with electrospray (voltage/flow rate) characteristics




References:  
Valaskovic, G.A, Murphy, J.P., Lee, M.S. "Automated Orthogonal Control System for Electrospray Ionization", J. Am. Soc. Mass Spectrom., 2004, 15, 1201-1215.



### Figure 11

Fourier transform analysis of MRFA 524 Da MH<sup>+</sup> ion at 2.0 kV spray voltage and 3.0 kV spray voltage. The plot of frequency (Hz) in the panels on the right indicate a much lower level of noise—indicating spray fluctuation—at the 3.0kV spray voltage. The decreased noise corresponds to more accurate peak area integration.



## Future Work

- Evaluate effect of spray stability on data quality (peak area) for quantitative calculations (RSDs)
- Investigate effect of spray stability on data quality at lower (200 nL/min) and higher (1-10 µl/min) flow rates
- Incorporate an orthogonal (frequency based) feedback control system via the digitally controlled nanospray source (Digital PicoView)